

Celestial Star™

The Celestial Star™ is ideal for surgical applications and procedures where maneuverability and compact design are key requirements. The light offers excellent performance in a compact package, very low maintenance and market leading reach and flexibility.

The Celestial Star™ mobile has significant advantages for mobile surgical applications. It's counter-balanced arm design makes it is extremely easy to move around.



CELESTIAL STAR™ FEATURES

- >>> Available in ceiling, dual ceiling, wall and mobile configurations
- >>> Low maintenance, space-saving design
- >>> 730 mm extension arm (optional 990 mm arm available) for maximum reach
- >>> 720° rotation allows head-to-toe coverage
- >>> Drift free balance arm with easy positioning
- >>> 3x35W long life Halogen bulbs (average 4,000 hours)



Mobile

TECHNICAL DATA

Power supply: 230-240V 50 Hz AC or 120V 60 Hz AC

Light intensity: Up to 64,500 Lux/6,000 Footcandles

Light source: Halogen bulb – Dichroic 12V/3x35W/10°

Light color temperature: Up to 4,100 Kelvin

Medical standards: FDA Class II

Color: White – RAL 9002

Configurations:	Ceiling	Dual Ceiling	Wall	Mobile
Weight:	23.6 kg	47,2 kg	17.3 kg	24.3 kg

Ceiling	Art. No.: CS2050C
Dual Ceiling	Art. No.: CS2050D
Wall	Art. No.: CS2050W
Mobile	Art. No.: CS2050M
Color Correction Filter	Art. No.: 18676

Note! 3 filters are required. Increases the color temperature to 4,100 Kelvin.



Wall



Duo Ceiling