

HM series – Multipurpose Examination lights

HM FEATURES

- >>> Available in wall, c-clamp, table and mobile configurations
- >>> Compact, vented shade minimizes heat output
- >>> Stay cool finger guard for easy adjustment of light
- >>> Swivel joint lamp head allows 359° axial rotation and 90° angular adjustment



TECHNICAL DATA

Power supply: 230-240V 50 Hz AC or 120V 60 Hz AC
Light intensity: Up to 27,000 Lux/2,500 Footcandles
Light source: Halogen bulb – 12V/20W
Light color temp: Up to 3,100 Kelvin
Color: White – RAL 9002

Configurations:	Ceiling	Wall	C-clamp	Table	Mobile	Mobile light	Floor
Weight:	2.8 kg	2.8 kg	0.9 kg	2.8 kg	7.2 kg	6.2 kg	3.8 kg

Ceiling	Art. No.: HM28007-2,5/2,7/3,0
Wall	Art. No.: HM28002W
C-clamp	Art. No.: HM741-20W
Table	Art. No.: 28013
Mobile, articulated arm	Art. No.: HM28009W
Mobile light, gooseneck arm	Art. No.: HM1006-37W
Floor	Art. No.: 28014

