

PF series – Multipurpose Examination lights

PF FEATURES

- >>> Available in ceiling, wall and mobile configurations
- >>> Swivel joint lamp head allows for easy vertical and horizontal adjustment of light
- >>> Adjustable arm levers allow for easy vertical and horizontal positioning
- >>> Choice of articulated arm or flexible gooseneck models



TECHNICAL DATA

Power supply: 230-240V 50 Hz AC or 120V 60 Hz AC
Light intensity: Up to 43,000 Lux/4,000 Footcandles
Light source: Halogen bulb – Dichroic 12V/50W/10°
Light color temp: Up to 3,000 Kelvin
Color: White – RAL 9002

Configurations:	Ceiling	Wall	Direct Mount	Mobile
Weight:	8.0 kg	7.7 kg	2.0 kg	10.4 kg

Ceiling	Art. No.: PF2017-44
Wall	Art. No.: PF2015-44
Direct Mount	Art. No.: PF28010
Mobile, articulated arm	Art. No.: PF2016-44
Mobile, gooseneck arm	Art. No.: PF2006-32

