

SoloSpot™

The SoloSpot™ is ideal for applications where a focused, intense light is required. The compact, specially designed lamp head offers cool light and a focused, sharp spot. This light is perfect for direct attachment to medical equipment such as patient tables and chairs.

SOLOSPOT™ FEATURES

- >>> Excellent maneuverability allows for easy adjustment of light
- >>> Intense 20W Halogen light source (average 2,000 hours)
- >>> 8,4 cm diameter spot pattern at 0,5 meters
- >>> Unique double filtered lens and color correction filter with lampshade extension helps keep light cool on skin and illuminate with natural skin tones and improve tissue clarity
- >>> Available in wall, rail, c-clamp and mobile configurations



TECHNICAL DATA

Power supply: 230-240V 50 Hz AC or 120V 60 Hz AC
Light intensity: Up to 25,000 Lux/2,300 Footcandles
Light source: Halogen bulb – 12V/20W
Light color temp: Up to 4,100 Kelvin
Color: White – RAL 9002

Configurations:	Wall	Rail	C-clamp	Mobile
Weight:	2.8 kg	2.8 kg	2.7 kg	7.1 kg

SoloSpot™ Art. No.: **28025**

Mounting Options

Wall	Art. No.: 26017
Rail	Art. No.: 28019
C-clamp	Art. No.: 18021
Mobile	Art. No.: 18468

