

TriStar™

Excellent performance and great value is what the TriStar™ has to offer. Ideal for medical procedures in environments where space is limited, and the range of procedures performed is broad. The unique telescoping pole offers height adjustment and the TriStar™ is one of the only products in the market with this functionality.

TRISTAR™ FEATURES

- >>> Available in wall and mobile configurations
- >>> Light pattern 125 mm diameter at 600 mm distance
- >>> Wall model offers 750 mm extension arm, providing 1400 mm reach to beam center
- >>> 3x20W long life Halogen bulbs (average 2,000 hours)
- >>> Mobile upright stand, raises from 685 mm to 2030 mm, via one hand control
- >>> Casters move easily across all surfaces; base design prevents tipping



Mobile

TECHNICAL DATA

Power supply: 230-240V 50 Hz AC or 120V 60 Hz AC
Light intensity: Up to 42,000 Lux/3,900 Footcandles
Light source: Halogen bulb – Dichroic 12V/3x20W/10°
Light color temperature: Up to 4,100 Kelvin
Medical standards: FDA Class II
Color: White – RAL 9002

Configurations:	Wall	Mobile
Weight:	8.2 kg	11.4 kg

Wall	Art. No.: TS2015W
Mobile	Art. No.: TS2015M
Color Correction Filter	Art. No.: 18676

Note! 3 filters are required. Increases the color temperature to 4,100 Kelvin.



Wall